

The Sexual Assault Specimen Collection.

Please contact the Microbiology Lab (x72585) for supplies and/or collection guidance.

Urine Cup	BD ESwab Collection and Transport System	APTIMA Multitest Swab Specimen Collection Kit
		

NOTE: Each collection site (except urines) should include the two swabs shown above. Please use both swabs for each respective site. Notify Microbiology of desired testing sites if there are any questions.

Test Names:				
	Chlam. trachomatis and N. gonorrhoeae PCR	Chlam. trachomatis and N. gonorrhoeae Aptima UW PCR	Trichomonas vaginalis PCR	
Specimen Source	Urine, Dirty Catch 1 Urine Cup	One Cup of Dirty Catch Urine		
	Throat	One BD E-swab	One APTIMA Multitest Swab	Not performed by UWMC
	Rectal	One BD E-swab	One APTIMA Multitest Swab	Not performed by UWMC
	Genital (Vaginal)	One BD E-swab	One APTIMA Multitest Swab	Not performed by UWMC
	Genital (Penile)	One APTIMA Multitest Swab (Sent to UWMC for testing)	One APTIMA Multitest Swab	Not performed by UWMC
	Other non-validated sites	One APTIMA Multitest Swab (Sent to Mayo Clinic for testing)	Not performed by UWMC	Not performed by UWMC

**** If clinically indicated for Herpes (HSV1 and 2), Test name: "HSV Rapid PCR, Swab" please call lab for further instructions.**

Lab use only test codes: CT GC PCR (LAB2394), GCCTAD (LAB2393), TRICAD (Order MISC Lab Test)

Last updated: 5/13/26