

The Sexual Assault Specimen Collection Kit. Please contact Microbiology Lab for pre-assembled kits (x72585)

Urine Cup	BD ESwab Collection and Transport System	APTIMA Multitest Swab Specimen Collection Kit
		

Each kit (except urines) will arrive with two swabs. Notify Microbiology of desired testing sites. Please use both swabs for each respective site.

Test Names:			
		Chlam. trachomatis and N. gonorrhoeae PCR	Chlam. trachomatis and N. gonorrhoeae Aptima UW PCR
Specimen Source	Urine, Dirty Catch 1 Urine Cup	One Cup of Dirty Catch Urine	
	Throat	One BD E-swab	One APTIMA Multitest Swab
	Rectal	One BD E-swab	One APTIMA Multitest Swab
	Genital (Vaginal)	One BD E-swab	One APTIMA Multitest Swab
	Genital (Penile)	One APTIMA Multitest Swab (Sent to UWMC for testing)	One APTIMA Multitest Swab
	Other non-validated sites	One APTIMA Multitest Swab (Sent to Mayo Clinic for testing)	Not performed by UWMC

**** If clinically indicated for Herpes (HSV1 and 2), Test name: "HSV Rapid PCR, Swab" please call lab for further instructions.**

Lab use only test codes: CT GC PCR (LAB2394), GCCTAD (LAB2393)

Last updated: 2/28/23