

The Sexual Assault Specimen Collection Kit. Please contact Microbiology Lab for pre-assembled kits (x72585).

Urine Cup	BD ESwab Collection and Transport System	APTIMA Multitest Swab Specimen Collection Kit	APTIMA Unisex Swab Specimen Collection Kit
			

Please collect the following swabs for STI testing according to specimen source:

		Tests		
		Culture, GC Screen + Wet Mount <u>C GCS+WM</u>	Chlam trachomatis and N. gonorrhoeae PCR <u>CT GC PCR</u>	Chlam trachomatis and GC Aptima UW PCR <u>GCCTAD</u>
Specimen Source	Urine, Dirty Catch 1 Urine Cup	one cup of Dirty Catch Urine		
	Throat 2 Types of Swabs	one BD ESwab		one APTIMA Multitest Swab
	Rectal 2 Types of Swabs	one BD ESwab		one APTIMA Multitest Swab
	Genital (Vaginal) 2 Types of Swabs	one BD ESwab		one APTIMA Multitest Swab
	Genital (Penile) 2 Types of Swabs	one BD ESwab		one APTIMA Multitest Swab
	Genital (Urethral) 2 Types of Swabs	one BD ESwab		one APTIMA Unisex Swab

If clinically indicated for Herpes (HSV1 and 2), call lab for COPAN Flocked Swab with UTM to collect for PCR Test (HSV Rapid PCR, Swab) from vaginal or penile sites.